



## 02MFSE00-1

## TNT SURGICAL MASK W AND VISOR SHIELD

- Disposable face mask made of four layers of non-woven fabric with anti-fog visor, adaptable nose clip
- With elastic headband, colour green
- Intended use: maintaining hygienic conditions between user-patient environment
- Made of TNT (non-woven fabric) is seamless
- The mask is non-toxic, hypoallergenic and odourless and offers very low breathing resistance
- It has a high filtering capacity and contains no glass fibres, nickel or other harmful substances
- · Latex free
- CE marking according to the European MDR 2017/745
- Reference standard and certification EN ISO 14683:2019 for type IIR masks
- Biocompatibility according to ISO 10993-1
- Expiry date: 3 years from the date of manufacture indicated on the packaging.
- Storage: in a cool, dry place away from strong heat sources

## Indications for use:

Wear the mask with clean, dry hands, with the side with the visor facing up and the support surface facing out. Then unroll the folds and shape the temple so that the mask fits as snugly as possible, making sure that the mouth and nose are adequately covered. Replace the mask if it is damaged, dirty or wet, or according to local regulations; do not touch the filter part of the mask after use and wash your hands after removing it.

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.