



**BH1511002MC**

## **STERILE SYRINGE 2.5 ML (100 PCS)**

- Sterile, single-use syringes with fitted needle and individually packaged
- Well-defined graduated scale. Body of high transparency.
- Synthetic rubber seal, free of latex components in order to exclude molecular transfer to the substances to be injected.
- Piston end stop to prevent the risk of spillage of the liquid to be injected or withdrawn .
- Free of latex and phthalates.
- Class IIa medical device in accordance with Directive 93/42/EC
- Size: 2.5 ml
- Pack of 100 pieces

*DIR/QM – Rev. 01\_11/08/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.