



**ENEBX2AY007**

## **RESCUE ISOTHERMIC BLANKET**

- Rescue and protection isothermal blanket
- It offers a dual function towards the patient's body, whether traumatized, burned or post-surgery:
  - By placing the silver part towards the body and the gold part on the outside, the blanket insulates and reflects while retaining about 80% of the body heat.
  - By placing the golden part towards the body and the silver part outside, the blanket reflects the sun rays and allows for the body to stay cool.
- Non-sterile
- Size: approximately 160x210 cm
- EC Medical device

*DIR/QM - Rev. 01\_26/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.