



**CER111**

## **PROTECTAPLAST SKIN-COLORED COHESIVE BANDAGE 3x450CM**

- Adhesive bandage to be used as a plaster on the wound
- It attaches to itself, without the need for adhesives
- Soft and elastic, easy to apply
- Allows the skin to breathe for better wound healing
- Absorbs blood and does not stick to the skin, hair or wound
- Can be used in combination with gel or cream or for fixing gauzes
- Easy to tear to size
- Latex-free
- Water-resistant
- Size: 3x450cm

*DIR/QM – Rev. 01\_31/03/2025*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.