



610004-20

## WRASFER - TRANSFER SHEET

- Highly absorbent transfer sheet
- Disposable, practical and hygienic, it reduces the costs of cleaning and sanitising stretchers
- Resistant: guarantees the transfer of patients up to 200 kg
- Able to absorb up to 5 litres of exudate: protects healthcare workers from contamination due to exudate leakage
- Sold in single non-sterile form and cannot be sterilised
- Intended for use in the health emergency sector, in hospitals, emergency rooms, nursing homes, mortuaries and for the recovery of corpses in special conditions of use. It can be used as a complement to other recovery and transport aids such as vacuum mattresses, vacuum stretchers, scoop stretchers and toboggan stretchers, etc.
- 3-layer structure: polyethylene, non-woven fabric and absorbent layer
- Maximum capacity: 200 kg
- Sheet dimensions: 220x100 cm
- Absorbent layer dimensions: 200x50 cm
- Weight: 350 g
- Colour: blue-white
- Risk class I medical device, compliant with Regulation (EU) 2017/745 (MDR).

**CAUTION:** the use of the device must involve at least two/four people, depending also on the patient's weight, and the conditions of use described in the device packaging must be respected.

The device cannot be used as a patient transport device (it is not a substitute for stretchers, beds or other medical devices intended for this purpose).

DIR/QM – Rev. 01\_03/04/2025

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.