



MIS. 3

2502466004

## TUBULAR ELASTIC NET SIZE 3 50M

- Elastic tubular net made of soft, stretchy yarn, suitable for securing dressings to medium-sized areas of the body.
- It optimally replaces traditional bandages
- Easy and quick application
- Dimensions: size 3 – 50 m tensed
- Size 3 indicated for: arm, leg, elbow, knee, shoulder

DIR/QM – Rev. 01\_06/03/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.