



MIS. 4

2502466005

## TUBULAR ELASTIC NET SIZE 4 50M

- Elastic tubular net made of soft, stretchable yarn, suitable for securing dressings on medium to large areas of the body.
- It optimally replaces traditional bandages
- Easy and quick application
- Dimensions: size 4 – 50 m tensed
- Size 4 indicated for: arm, leg, elbow, knee, shoulder, head, thigh, hip

DIR/QM – Rev. 01\_06/03/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.