



811023E

NEEDLE FIXATION DRESSING 6X8 CM

- Self-adhesive dressing for stabilising, protecting and securing intravenous cannulas
- Features:
 - A recess for vertical insertion of the intravenous cannula
 - A central absorbent pad that prevents the non-woven fabric from sticking to the injection site
 - Rounded edges that prevent accidental removal of the dressing
 - Silicone-coated paper backing that facilitates precise, painless and efficient application
 - Made from hydrophobic non-woven fabric
 - Hypoallergenic
 - Shelf life: 5 years
 - Size: 6 cm x 8 cm
 - Pack of 50

DIR/QM – Rev. 01_08/06/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.