



FP6100

FABRIC PLASTER 6CM x 1M

- Device intended for the protection of minor wounds, skin lesions, cuts, grazes and lacerations with mild to moderate exudate.
- Made from cotton fabric with a central absorbent pad, capable of absorbing small to medium amounts of exudate and protecting the wound from external agents.
- Flexible and conforms to the contours of the body, ensuring comfort and freedom of movement.
- Equipped with a silicone-coated paper backing to allow for precise, simple and non-traumatic application.
- Good permeability to air and water vapour.
- Skin-coloured.
- Single-use.
- Non-sterile.
- Dimensions: 6 cm x 1 m.
- Individually packaged.

DIR/QM – Rev. 01_28/05/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.