



ST-002

TACTICAL TOURNIQUET

- Made from soft, breathable material to prevent chafing and skin burns
- Features:
 - A dual-closure system that eliminates the need for screws and clips, and a locking mechanism to prevent the tourniquet from opening accidentally
 - One-handed operation
 - Durable loop system with a movable inner strap that provides circumferential compression to the limb
 - Additional element to secure the tourniquet and control blood flow
 - Space to note the time of application
- Individually packaged
- Latex-free, phthalate-free
- Disposable and non-sterile
- Weight: 82 g
- Dimensions when closed: 16 cm x 6 cm x 4 cm
- Dimensions when open: 96 cm x 4 cm

DIR/QM – Rev. 01_20/05/2026

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.