



**BOR84064**

### **E-HOLSTER BLUE COMPLETE\*\*\***

- Complete belt E-Holster
- Made in tearproof and water-resistant blue material
- 5 compartments with pockets for gloves and for Pen Light
- Belt loop of about 5,5 cm
- Dimensions: about 14x3x19(h) cm
- Complete with:
  - 1 first aid scissors
  - 1 anatomic forceps
  - Cocher forceps
  - 1 Pen Light
  - 1 Tourniquet

*DIR/QM – Rev. 02\_08/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.