

## **Emergency Medical Systems**







## CCPE1200PER

## PERFECTO PE SHOE COVER 100PCS

- Blue disposable shoe cover in PE (polyethylene) 10 μm.
- Widely used in hospitals, clinics, beauty salons, food industries to avoid carrying dirt from the outside to the inside.
- Medical Device Class I according to European Directive 93/42/EEC and following amendments (including Directive 2007/47/EC).
- Lot number indicated on the box.
- Storage: Room temperature in a dry place, away from heat sources, flames, or sparks.
- Packaging: Box of 100 pcs.

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.