



**ACC180018**

## **LIGHT BLUE NURSE WATCH**

- Silicone watch for nurses and medical personnel
- Equipped with a pin to attach it to the uniform
- Suitable for measuring patient pulse rates and times for sterilization procedures
- Complies with current regulations prohibiting nurses from wearing wristwatches
- Single-color round model
- Dimensions: watch face 35 mm, watch length 85 mm
- Weight: about 30 g
- Battery: SR626SW button cell (battery life 18 months)
- Not a medical device

*DIR/QM – Rev. 01\_22/05/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.