



**ACC180018**

## **LIGHT BLUE NURSE WATCH**

- Silicone watch for nurses and medical personnel.
- Equipped with a safety pin for attaching to the healthcare uniform.
- Suitable for measuring patient pulse rates and times for sterilization procedures.
- Complies with current regulations prohibiting nurses from wearing wristwatches.
- Single-color round model.
- Dimensions: watch face 35 mm, watch length 85 mm.
- Weight: about 30 g.
- Battery: SR626SW button cell (battery life 18 months).
- Not a medical device.

*DIR/QM – Rev. 01\_22/05/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.