



## 16100

## **BASKET RESCUE STRETCHER**

- Rescue stretcher, extremely lightweight, suitable for patient transport in uncomfortable environment conditions
- Made in polymeric ABS material with aluminium frame
- Complete of expanded polyutherane mattress and adjustable footrest
- 4 quick unhook security belts
- 12 transport handles
- Arranged for the anchorage to helicopters
- Dimensions: about 216x61x17h cm
- Weight: about 15 kg
- Capacity: max 300 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.