



9031

## CHAMPION COLLAR - PEDIATRIC

- Emergency Cervical Collar Champion Collar is new double valve collar adjustable in different heights.
- Allows immobilization of the cervical spine, minimizes flexion-extension and ensures better control of axial rotation and lateral flexion movements.
- Ideal for those who want to have, with a single article, several possibilities of immobilization.
- Ideal for post-traumatic situations in the cervical area, where a valid support and/or protection of the skeletal system and the muscle bands is necessary.
- Central hole for the control of the trachea and the execution of a tracheotomy procedure.
- Has been developed for use during medical emergency and flat two-piece design allows efficient placement in limited spaces.
- The extremity of the Velcro is in orange color for a better identification of the grip.
- Polyethylene structure, black LD padding, nylon buttons and Velcro strap latch.
- If used when skin lesions occurred, protect the wound with an appropriate bandage; in such situations, do not use the device in direct contact with the wound.
- Can be worn during CAT scan and MRI.
- 100% Radiotransparent - Latex-Free - Unsterile - Single-use ? Disposable
- The device is disposable. However, it is possible to wash in water at 30°C with neutral soap or detergent.
- RDM Code 1424420/R
- CND Code Y060312

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.