



BOR84064

E-HOLSTER BLUE COMPLETE***

- Complete belt E-Holster
- Made in tearproof and water-resistant blue material
- 5 compartments with pockets for gloves and for Pen Light
- Belt loop of about 5,5 cm
- Dimensions: about 14x3x19(h) cm
- Complete with:
 - 1 first aid scissors
 - 1 anatomic forceps
 - Cocher forceps
 - 1 Pen Light
 - 1 Tourniquet

DIR/QM – Rev. 02_08/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.