



**IMM120060**

## **BABY VACUUM MATTRESS**

- BABY vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of children
- Equipped with:
- 4 handles for lifting and carrying
- Suction valve located close to the head on the outer side, easily accessible even with immobilized patient
- Radio-transparent
- Dimensions: 87x65x4h cm
- Weight: 1.6 kg
- Max. load capacity: 18 kg
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

*DIR/QM – Rev. 02\_21/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.