















## IMM120071

## **INFANT TRANSPORT UNIT**

- Newborn emergency transportation system
- Complete with inner fixing belts
- · Grants a warm, protected and safe environment
- Innovative foldable-bag system cordura material allowing, when not used, to be taken everywhere, always ready for use
- 3 transparent sides allowing complete and continuous monitoring of newborn by medical personnel
- Internal mattress ensuring comfortable transport
- 2 external handles for transportation
- Arranged for the anchorage to the ambulance stretcher (EN 1789:1999)
- Dimensions open: cm 90x38x34,5(h)
- Dimensions closed: cm 90x38x15(h)
- Supplied with transport bag

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.