



**IMM120129**

## **PROFI EVOLUTION VACUUM MATTRESS-**

- PROFI EVOLUTION vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients
- Equipped with:
- semi-rigid support in which 8 handles for lifting and transport are located
- suction valve located close to the head on the outer side, protected by the semi-rigid support, easily accessible even with immobilized patient
- 2 specific holes for fixing the mattress to metal stretchers used in air rescue operations
- 3 immobilisation straps
- 3 internal suction chambers for better adaptation to the anatomical shape of the patient
- Radio-transparent
- In compliance with the standard EN 1865-1
- Dimensions: 205x95x10h cm
- Weight: 8.2 kg
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

DIR/QM – Rev. 02\_21/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.