



IMM120150

PROFI VACUUM MATTRESS

- PROFI vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients
- Equipped with:
 - semi-rigid support in which 8 handles for lifting and transport are located
 - suction valve located close to the head on the outer side, easily accessible even with immobilized patient
 - 3 immobilisation straps
 - 3 internal suction chambers for better adaptation to the anatomical shape of the patient
 - Radio-transparent
 - In compliance with the standard EN 1865-1
 - Dimensions: 205x95x10h cm
 - Weight: 10 kg
 - Maximum load capacity: 150 kg
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

DIR/QM – Rev. 02_21/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.