



**IMM120150**

## **PROFI VACUUM MATTRESS**

- PROFI vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients
- Equipped with:
  - semi-rigid support in which 8 handles for lifting and transport are located
  - suction valve located close to the head on the outer side, easily accessible even with immobilized patient
  - 3 immobilisation straps
  - 3 internal suction chambers for better adaptation to the anatomical shape of the patient
  - Radio-transparent
- In compliance with the standard EN 1865-1
- Dimensions: 205x95x10h cm
- Weight: 10 kg
- Maximum load capacity: 150 kg
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

*DIR/QM – Rev. 02\_21/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.