



## IMM120172

## **VACUUM SPLINTS FOR LEG 108CM**

- Made of high-strength, easily washable fabric for quick and painless immobilisation of all types of fractures, even compound fractures
- Immobilises but does not compress and allows natural blood circulation
- Can also be used on children
- Equipped with a plastic suction valve on the outer side
- Radiotranslucent
- Dimensions: approx. 108x64 cm
- Weight: approx. 1.7 kg
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the eIFU webpage

DIR/QM - Rev. 02\_21/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.