



**IMM120200**

## **2 VACUUM SPLINTS SET + BAG**

- Set of 3 vacuum splints: leg and arm
- Suitable for paediatric use
- Absolute, fast and painless immobilisation of fractured limbs (even exposed ones)
- High resistant washable material
- Compact suction plastic valve with external  $\frac{1}{4}$  on/off.
- Ensuring a perfect immobilisation but do not compress, and allow physiological blood circulation
- 100% X-ray transparent (no metal parts included)
- Complete with transport bag and mini-pump
- Arm splint dimensions: about 47x75 cm
- Leg splint dimensions: about 64x92 cm
- Dimensioni of set: about 54x20x44(h) cm
- Weight: about 3,2 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.