



## IMM121640

## **UNIVERSAL HEAD IMMOBILIZER - WIND**

- Realized in stuffed MDF and covered with washable, disinfectable, water repellent and thermo-welded material
- 2 head-fixing-pillows with opening for access to the ears are integrated at the base and are part of the base itself if not used
- Allows fixed head immobilizer storage with stretcher
- Complete of 2 chin straps and 2 Velcro straps, 8 belts with loop and 2 velcro bands for the fixing to the spinal board.
- X-ray translucent.
- Dimensions: about 44,5x30x2,5 cm
- Weight: about 1,2 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.