



IMM124034

SPLINTART WRIST SSPLINT

- Splint with aluminium core for easy limbs immobilization in case of fracture or sprain
- Made in neoprene material
- Adaptable to fit each body frame
- It does not produce tissues compressions
- It maintains body temperature
- It reduces congestion around injured areas
- Velcro modelling closing
- Lightweight and easily washable even in washing-machine
- Least encumbrance
- Dimensions: about 32x15x0,6 cm
- Weight: about 150 gm
- Not suitable for femur fractures

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.