



**IMM124035**

## **ELBOW/ANKLE SPLINTART SSPLINT**

- Splint with aluminium core for easy limbs immobilization in case of fracture or sprain
- Made in neoprene material
- Adaptable to fit each body frame
- It does not produce tissues compressions
- It maintains body temperature
- It reduces congestion around injured areas
- Velcro modelling closing
- Lightweight and easily washable even in washing-machine
- Least encumbrance
- Dimensions: about 57x18x0,6 cm
- Weight: about 280 gm
- Not suitable for femur fractures

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.