



**IMM121033**

## **FIXING SYSTEM FOR SPINAL STRETCHERS**

- Allows the immobilization of an adult or paediatric patient on a generic spine board
- The belts are made with 1.2 mm thick polyester webbing to ensure high fastening strength and durability
- Thanks to the possibility of moving the straps according to the type of patient, (along the vertical axis of the fastening system) it is possible to secure the patient in a stable and safe way
- The straps are fastened to the patient's body thanks to 5 cm wide Velcro® strips (approximately) which make both closing and opening operations easy and safe
- The system is complete with a practical case
- The device can only be washed with cold water. Do not use aggressive solvents or detergents.

*DIR/QM – Rev. 02\_26/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.