



**LD-588**

### **AUTOMATIC/DIGITAL SPHYGMOMANOMETER\***

- 52x59 mm LCD display
- Up to 60 measurements storage capacity
- Oscillometric measurement method
- Pressure measurement scale: 40-260 mm/Hg
- Pulse measurement scale: 40-160 beats/min
- Accuracy:  $\pm 3$ mm/Hg (pressure)  $\pm 5\%$  (pulse)
- Automatic pressurization
- Quick depressurization with valve
- Power supply: 4 x 1.5V AA batteries (not included)
- Dimensions: 127x99x50 mm
- Weight: approx. 255 g without batteries
- It is a CE medical device
- Measurements are to be carried out by qualified healthcare professionals

*DIR/QM – Rev. 02\_04/09/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.