



## KED PED

### PAEDIATRIC SPINAL IMMOBILIZER

- Complete with pillow and carrying bag
- Allows the immobilization of trunk, head and spine and extraction of patients from difficult to access areas
- Wrap-around and flexible design
- Equipped with coloured ribbon and quick release/hook buckles
- Equipped with shaped wooden splints that give rigidity to the device
- Made of washable blue fabric
- Radiolucent
- Class I medical device in compliance with Regulation (EU) 2017/745
- Size: paediatric

DIR/QM – Rev. 02\_28/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.