



XIMM120014

SDM VACUUM MATTRESS COMPLETE

- SDM vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients.
- Equipped with:
- 8 handles for lifting and transportation
- suction valve located close to the head on the outer side, easily accessible even with immobilized patient
- 3 immobilisation straps
- 3 internal suction chambers for better adaptation to the anatomical shape of the patient
- Radio-transparent
- Complete with suction pump and carrying bag
- Dimensions: 200x80x6h cm
- Weight: 7 kg
- In compliance with the standard EN 1865-1
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the eIFU webpage

DIR/QM - Rev. 02_28/08/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.