



**XIMM120014**

## **SDM VACUUM MATTRESS COMPLETE**

- SDM vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients.
- Equipped with:
  - 8 handles for lifting and transportation
  - suction valve located close to the head on the outer side, easily accessible even with immobilized patient
  - 3 immobilisation straps
  - 3 internal suction chambers for better adaptation to the anatomical shape of the patient
  - Radio-transparent
  - Complete with suction pump and carrying bag
  - Dimensions: 200x80x6h cm
  - Weight: 7 kg
  - In compliance with the standard EN 1865-1
  - It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

*DIR/QM – Rev. 02\_28/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.