



**XIMM120130**

## **PROFI VACUUM MATTRESS COMPLETE**

- PROFIL vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients
- Equipped with:
  - a semi-rigid support in which 8 handles for lifting and transport are located
  - Suction valve located close to the head on the outer side, easily accessible even with immobilized patient
  - 3 immobilisation straps
  - internal suction chambers for better adaptation to the anatomical shape of the patient
  - Radio-transparent
- Complete with suction pump and transportation bag
- Dimensions: 205x95x8h cm
- Weight: 10 kg ca.
- In compliance with the standard EN 1865-1
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

*DIR/QM – Rev. 02\_28/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.