



XIMM120130

PROFI VACUUM MATTRESS COMPLETE

- PROFI vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients
- The mattress includes a semi-rigid support in which 8 handles for lifting and transport are located
- Suction valve located close to the head on the outer side, easily accessible even with immobilized patient
- Equipped with 3 immobilisation straps
- Radio-transparent
- Equipped with 3 internal suction chambers for better adaptation to the anatomical shape of the patient
- Equipped with suction pump and transportation bag
- Dimensions: 205x95x8h cm
- Weight: 10 kg ca.
- In compliance with the standard EN 1865-1
- It is a CE medical device in compliance with the European Regulation MDR 2017/745

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.