



**XIMM120131**

## **PROFI VACUUM MATTRESS EVOLUTION COMPLETE**

- PROFI EVOLUTION vacuum mattress, made of high-strength fabric, intended for the immobilisation and transport of polytrauma patients
- Equipped with:
  - A semi-rigid support in which 8 handles for lifting and transport are located
  - Suction valve located close to the head on the outer side, protected by the semi-rigid support, easily accessible even with immobilized patient
  - 2 specific holes for fixing the mattress to metal stretchers used in air rescue operations
  - 3 immobilisation straps
  - 3 internal suction chambers for better adaptation to the anatomical shape of the patient
  - Radio-transparent
  - Complete with suction pump and transportation bag
  - Dimensions: 205x95x10h cm
  - Weight: 8.2 kg
  - In compliance with the standard EN 1865-1
  - It is a CE medical device in compliance with the European Regulation MDR 2017/745

User manual available on the [eIFU](#) webpage

*DIR/QM – Rev. 02\_28/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.