



**IMM120214**

## **ARM VACUUM SPLINT BASIC+TRANSPORT BAG**

- Anatomic conformed vacuum splint suitable for leg and ankle
- Fast and painless immobilisation of all types of fractures, even compound fractures
- Immobilises but does not compress and allows normal blood circulation
- High-strength fabric, easily washable
- Equipped with:
  - plastic suction valve positioned on the outer side
  - Tear-proof and waterproof fabric protection bag with quick-release fastening system
  - Radio-transparent
- Dimensions: approx. 68x38x53 cm
- Weight: approx. 0.8 kg
- EC Medical device

*DIR/QM – Rev. 02\_21/08/2024*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.