



**IMM121631R**

## **UNIVERSAL OB HEAD IMMOBILIZER RED**

- Made of treated, water-repellent high-density foam in red color
- Consists of 1 base cushion, 2 head support pads with opening for ear access, 2 chin rests with central padding and 8 straps, 2 tape straps with Velcro already integrated at the base for positioning the head immobilizer on the stretcher
- Radio-transparent
- Dimensions: 40x26x16 cm approx.
- Weight: 880 g approx.

*DIR/QM – Rev. 01\_15/12/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.