



**IMM121631P**

## **UNIVERSAL OB HEAD IMMOBILIZER PEDIATRIC**

- Made of treated, water-repellent high-density foam in red colour
- Consisting of 1 base cushion, 2 head support pads with opening for access to ears, 2 chin rests with central padding and 4 straps, 3 Velcro straps already integrated in the base for positioning the head immobilizer on the stretcher
- Radio-transparent
- Dimensions: approx. 34x26x15 cm
- Weight: approx. 500 g

*DIR/QM – Rev. 01\_15/12/2023*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.