



YHR-ED01

EXTRICATION DEVICE

- Complete with pillow and carrying bag
- Allows immobilisation of trunk, head and spine and extraction of patients from hard-to-access areas
- Wrap-around and flexible design
- Equipped with coloured ribbon and quick release/hook buckles
- Equipped with shaped splints that give rigidity to the device
- Made of washable green fabric
- Radiolucent
- Class I medical device in accordance with Regulation (EU) 2017/745
- Size: 85x95x57 adult

DIR/QM – Rev. 01_31/03/2025

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.