



**YHR-CS1**

## **SCOOP STRETCHER**

- Aluminium scoop stretcher
- Effective hook/unhook device for stretcher parting
- Adjustable length with automatic block
- Bag with n. 3 quick unhook safety belts for patient immobilization
- Capacity: about 200 Kg
- Weight: about 8,5 Kg
- Maximum Size (LxWxH)?210x44x6cm
- Minimum Size (LxWxH)?120x44x9cm

*DIR/QM – Rev. 01\_31/03/2025*

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.