



BOR87002N

EMERGENCY RUCKSACK MINI VIPER 2.0

- Rescue Rucksack Viper Mini 2.0
- Red tearproof and water-proof material
- Equipped with double upper and side handle for easier transport
- Large front opening pocket with 2 internal elastic bands for the fixing of material
- 2 small side pockets and a large one on the outside
- Equipped with front, side and rear reflective strips (on padded shoulder straps)
- The part in contact with the back of the user is equipped with cushions
- 4 removable coloured pouches with transparent front and Velcro strips
- Reinforced bottom and rubber feet
- Dimensions: about 44x26x44cm
- Weight: about 2.3kg
- Maximum weight capacity: 30kg
- It can be washed by hand at a maximum temperature of 30° C
- Supplied empty but can be equipped with medical devices upon customer request
- Rescue Rucksack Viper Mini complies with REACH regulations but is not a medical device

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.