



**BOR89010**

## **PROFI BAG EMPTY**

- Emergency empty Profi Bag
- Made of impermeable and tearproof red fabric material
- With handle and shoulder straps
- Modular interior
- 1 frontal pocket
- Arranged to hold 1-2 lts oxygen cylinder
- 2 wide lateral external pockets
- Internal pocket with transparent.
- Dimensions: about 60x25x27(h) cm
- Weight: 1,4 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.