



**BOR89050**

### **TRAUMA BAG 20 RED EMPTY\***

- Functional, capacious, resistant and compact emergency bag
- Made of tearproof and impermeable red fabric material
- 4 refractive straps
- Totally removable cover
- Easy to carry either by hand, on shoulder or as a backpack
- 2 external lateral pockets
- 2 external lateral pockets with elastics and additional pockets
- 5 inner modular spaces
- 1 inner removable pouch with transparent frontal
- Dimensions: about 53x33x34(h) cm
- Weight: about 2,2 kg
- 

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.