



DIA01012Y

STETOSCOPIO SINGLE HEAD GIALLO*

- Compact and lightweight blue flat-head stethoscope
- Yellow “Y” tube made of a single piece (PVC material)
- Very sensitive membrane with excellent tone transmission
- Intended for all sectors of medicine (emergency, first aid, clinic, and hospital)
- CE marked according to the recent legal provisions concerning medical devices
- For non-professional and professional use
- Risk class according to Directive 93/42/CEE: I
- The device complies with the REACH regulation in force

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.