



**DIA02009**

### **ANEROID TABLE SPHYGMOMANOMETER\***

- Circular wall and/or table aneroid sphygmomanometer (diameter 170mm)
- Scale from 0 to 300 mm/Hg
- Measuring accuracy:  $\pm 3$  mmHg
- Velcro-closing cuff and inflation bulb, spiral connecting tube
- Wall anchor support with basket bracelet holder
- For professional and non-professional use
- Risk class I m, according to directive MDD 93/42/CEE (measuring instrument)
- The device complies with the REACH regulation in force

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.