



HS10A

HANDHELD PULSEOXIMETER

- Lightweight for carrying and Easy-To-Use
- Silicon rubber shell protection and stable metal bracket for table usage
- Adult, pediatric and neonate oximeter probe available
- Class II protection according to the IEC60601-1
- Protection against electric shock: type BF applied part
- Protection against the penetration of fluids and substances: IPX1
- Visual and acoustic alarm available according to the IEC60601-1-8
- 45 to 85 dB adjustable acoustic alarms
- Big size 2.4 LCD display with SPO₂, pulse ratio, pulse bar-graph and plethysmogram
- Friendly easy to set menu to adjust all parameters
- SPO₂ resolution 1 %
- Low Battery voltage indicator
- Refreshing rate: 1 sec
- The unit can work with alkaline batteries (2000 mAh) or rechargeable Lion battery pack (1600 mAh)
- Run time with new alkaline batteries: 30 hours without acoustic alarms and brightness set to the minimum at the temperature of 25° C
- Run time with fully charged Lion battery: 24 hours without acoustic alarms and brightness set to the minimum at the temperature of 25° C
- Charger stand available upon request
- User manual available in English, Italian and German languages
- Warranty: 2 years from the purchasing date

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.