



#### MED98874

### REINTEGRATION FIRST-AID BOX COMPL.

Medication package designed for companies or manufacturing units of group A and B and in any case with three or more workers. Content according to DM 388/2003 attachment 1 and D.L. 81/2008 art.45. For group A companies, the content provided for in attachment 1 is to be integrated on the basis of the risks present in the workplace.

complete with:

- 5 pairs of disposable sterile gloves;
- 1 anti-squirt visor;
- 1 litre of cutaneous Povidone-iodine with 10% iodine for professional use;
- 1.5 litre of saline solution (sodium chloride – 0,9%);
- 10 packages of sterile gauzes 10 x 10 cm – single packaging;
- 2 packages of sterile gauzes 18 x 40 cm – single packaging;
- 2 disposable sterile sheets;
- 2 disposable sterile medical forceps;
- 1 medium-sized elastic bandage;
- 1 package of hydrophilic cotton;
- 2 packages of various sized plasters ready-to-use;
- 2 rolls of high plaster 2,5cm;
- 1 pair scissors;
- 3 tourniquets;
- 2 packages of instant ice;
- 2 disposable sanitary waste bags;
- 1 thermometer;
- 1 sphygmomanometer;
- 1 stethoscope;
- Multilingual First Aid instructions;
- Copy of DM (ministerial order) 388/2003
- Dimensions: about 40x26x13 (h) cm
- Weight: about 4.6 kg

(Group A companies are those at high risk, such as thermoelectric power plants, nuclear plants and laboratories, mining

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.



companies, explosives factories, or companies/production units with more than 5 workers with high INAIL accident index or companies/production units with more than 5 permanent workers in the agricultural sector).

(Group B companies are those with 3 or more workers that do not fall into Group A, so they are not high risk).

(Group C companies are those with less than 3 workers that do not fall into group A, so they are not at high risk).

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.