



VAL44410

PAEDIATRIC INTUBATION SET IN POUCH*

- Complete intubation set in pouch with transparent frontal, featuring:
- 1 laryngoscope pediatric handle
- 1 C.L. Mc Intosh blade size 0
- 1 C.L. Miller blade size 1
- 1 C.L. Miller blade size 0
- 1 Magill forceps, size 16
- 1 Pean forceps
- 4 endotracheal tubes
- 4 Guedel airways sizes 000/00/0/1
- 1 intubation stylet
- 1 plaster roll
- Dimensions: about 27x12x9(h) cm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.