



**VBM30-04-0041**

## **ADULT QUICKTRACH**

- Urgent CONIOTOMY sterile set QUICKTRACH type
- Complete with one plastic canula with intubation internal-blade stylet with diameter of 4 mm
- Complete with ISO joint for the connection to resuscitation bag
- Quicktrach features a sharp conical-tip intubation stylet that avoids windpipe incision using classic scalpel and a ?lock? that avoids to penetrate too much in depth
- 2 lateral wings allowing airways safe fixing

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.