



#### OXIP7004

### OXI CYLINDER SWITCH 2IN-1OUT

- Mechanical manual switcher for oxygen control on oxygen distribution kits with 2 inlets and 1 output
- Both available with lateral and back inlets/outlets
- It allows 2 low pressure oxygen flow selections (max. 8 bar – 800 kPa)
- Oxygen compatible, made with properly treated material on surface so to grant full compliance to ISO15001
- Class of risk: IIa
- Big control knob easy to use also with gloves
- 3/8" oxygen hose connector
- Optional: frontal plate for wall embedding

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.